

LD-500

Automated HbA1c Analyzer



Sample Type	Whole blood and Prediluted	Interface	RS-232C, USB
Parameters	HbA1c	Operating Temperature	15-30°C
Principles	HPLC	Display	7" color touch screen
Wavelength	415nm	Printer	80mm integrated thermal
Time Per Test	< 5 minute	Autoloader Capacity	40tubes, continuously loading available
Testing Range	3-18%	Power Supply	AC 110/220V, 50/60Hz
Precision	CV<=3%	Dimension	W610mm×L435mm×H420mm
Sample Volume	10ul	Weight	Approximately 24kg

Accomplish **More**
In Less Steps



LABNOVATION TECHNOLOGIES, INC.

Address: 1F, W2-B Building, Keji 2 Road South, Hi-Tech Industrial Park Shenzhen 518057, China
Phone: +86 755 8636 8328 and +86 755 2652 5802
Fax: +86 755 8636 8318
Email: export@labnovation.com
Web: www.labnovation.com

